

AD-A268 802



DATA ACQUISITION FROM BLAST  
OVERPRESSURE TRIALS

AR-008-234


D.R. KIRK

MRL-TN-628

MARCH 1993

DTIC  
ELECTE  
AUG 30 1993  
S E D

APPROVED  
FOR PUBLIC RELEASE



MATERIALS RESEARCH LABORATORY

DSTO

# *Data Acquisition from Blast Overpressure Trials*

*D.R. Kirk*

MRL Technical Note  
MRL-TN-628

## *Abstract*

*A Macintosh computer has been used to acquire data from blast overpressure trials on various weapons. The computer is connected to a multiple channel FM data recorder via a MacSCSI488 bus controller, allowing the computer to control the recorder and to acquire data from it through an analog to digital converter. Detailed instructions are given for connecting the hardware and operating the software involved.*

DSTO MATERIALS RESEARCH LABORATORY

93 8 26 07 9

93-20051  
■■■■■■■■■■

*Published by*

*Materials Research Laboratory  
Cordite Avenue, Maribyrnong  
Victoria, 3032 Australia*

*Telephone: (03) 246 8111*

*Fax: (03) 246 8999*

*© Commonwealth of Australia 1993*

*AR No. 008-234*

**APPROVED FOR PUBLIC RELEASE**

## Contents

1. INTRODUCTION 5
  2. SOFTWARE 6
  3. HARDWARE CONNECTIONS 6
  4. STARTING UP 7
  5. THE DATA ACQUISITION PROGRAM 8
  6. OPERATING THE ACQUISITION PROGRAM 9
  7. MANIPULATING THE DATA 11
  8. CONCLUSION 14
- APPENDIX I — *Listing of Program for One Channel Acquisition* 15
- APPENDIX II — *Listing of Program for Two Channel Acquisition* 18

Accession For	
NTIS CRA&I	<input checked="checked" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification .....	
By .....	
Distribution /	
Availability Codes	
Dist	Avail and/or Special
A-1	

DTIC QUALITY INSPECTED 8

# *Data Acquisition from Blast Overpressure Trials*

## *1. Introduction*

The firing of guns and other explosive weapons produces a rapid change in pressure in the atmosphere. These pressure changes can damage the hearing of any person close to the weapon, particularly the operator. To assess levels of overpressure, trials are performed at various sites, depending on the nature of the weapon being tested. The pressures are measured at a number of locations using pressure transducers to convert pressure to voltage which is recorded on a high speed FM data tape recorder. A typical recorder can store 14 channels of data at a speed of 76 cm/s. A playback speed of 1.2 cm/s gives a time magnification of about 63. Thus 1 ms of real time expands to 63 ms of tape time, allowing for much greater detail to be examined.

Records from the FM data recorder can be transferred to an oscilloscope and plotted with a pen recorder. Alternatively, the tape can be played through an analog to digital (A-to-D) converter and data transferred to a computer for a more detailed analysis. An arrangement that has been used for some time in Explosives Ordnance Division, Salisbury consists of transferring the data to a Macintosh computer via a MacSCSI488 Iotech Bus Controller. This is simply a device that converts IEEE488 interface signals to Small Computer Systems Interface (SCSI) signals and vice versa, thus allowing the Macintosh to communicate with IEEE devices, such as a data recorder and an A-to-D converter. The advantage of the Macintosh is the variety of software that is available to manipulate data and the convenient user interface. The software suite currently used to produce final pressure-time data includes a locally written program using Microsoft Quickbasic, to effect the transfer of data and also to manipulate the data acquired; McSink, a Macintosh desk accessory that allows hidden characters (such as carriage controls) to be edited; and Igor, a data manipulation and graph plotting program from WaveMetrics.

---

## 2. Software

A program written in Microsoft QuickBasic is used to transfer the data from the FM data recorder to the Macintosh via a 12-bit analog to digital converter. The program takes into account the voltage range used to record the data, the pressure to voltage conversion factor and the record and playback speeds. It prompts the user for information such as the file name to store a channel of data. It controls the operation of the recorder and the A-to-D converter, such as start and stop. All the operator has to do is to position the tape at an appropriate starting point and activate the program.

The data acquired from the FM data recorder are in a raw form and consist of columns of figures for time and pressure and a third column if a trigger channel is recorded. The zero point for the time base is arbitrary and if a time of arrival of a certain data value, for example maximum pressure, from a given point is required the data must be manipulated to reflect this. The WaveMetrics program "Igor" is ideal for this purpose. Data can be displayed and edited in tabular and graphic form. Factors can be applied to the data to adjust the time base to the correct zero point and a preliminary plot can be displayed on the screen. A Macintosh Desk Accessory McSink is used at a later time to further transform the data.

Some software is also required to drive the MacSCSI488 Interface. This is supplied by Iotech and must be installed according to the product manual.

## 3. Hardware Connections

Several pieces of hardware are needed to process the data. A Kyowa 14 channel RTP-650A data recorder is commonly used but any recorder can be substituted provided it has a minimum of 40kHz bandwidth, FM mode, multiple speeds and IEEE (GPIB) Interface. An analog-to-digital converter such as the Kyowa ADC-116C is used to convert  $\pm 5$  volts signed input to a digital range of - 2048 to + 2048, an 11 bit range, with the twelfth bit indicating the sign. An Iotech MacSCSI488 bus controller converts IEEE signals to SCSI signals and vice versa to allow a Macintosh computer, typically an SE/30, to communicate with and control the recorder and the A-to-D converter. Subsequent discussion in this report applies to the hardware described but adaptations can be made for other hardware. Two IEEE488 cables are needed to connect the MacSCSI488 to the A-to-D converter and the recorder and a SCSI cable is used to connect the Macintosh to the MacSCSI488 bus controller. The data are transferred through one or two coaxial cables. To avoid the continual connecting and disconnecting of cables the pressure data are usually acquired by connecting the MONI output on the recorder to channel 1 of the A-to-D converter. This allows the operator to select the channel to be monitored by using the button on the recorder labelled "CHANNEL SELECT". If a trigger channel was recorded this channel is connected to the channel 2 input of the A-to-D converter. A diagram of the cable connections required is shown in Figure 1. Sometimes the trigger peak is a very narrow pulse and does not always register on the A-to-D converter, especially if the sampling rate is relatively slow (200 Hz for example). A small locally built electronic circuit is used to detect this peak and prolong it sufficiently to be recognised by the A

to D converter. This device is connected between the recorder and the A to D converter in the trigger channel coaxial cable.




Figure 1: Cable connections for transfer of data from VTR to Macintosh.

#### 4. Starting Up

Before acquiring any data it has been found useful to record a copy of the original tape and use that in the recorder, rather than the original. A faulty recorder may damage the tape and data may be lost. A copy can be made by simply connecting the channel outputs from one data recorder to the channel inputs of another data recorder and playing the original tape on one recorder while recording on the other.

When starting up the various devices it is important to switch on the data recorder and the A-to-D converter first, then the MacSCSI488 bus controller, then the computer. The mouse button should be held down when switching on the computer and held until a list of INITs appears. One is called MacSCSI488 which must be selected as active. There must be no external drives on the SCSI interface, and any associated INITs should be switched off. It is also a good idea to disable any screen savers, which can cause loss of data on screen when the computer is acquiring data. Other INITs which might be disabled include TOPS and TOPS Spool but virus protection software such as Gatekeeper and Disinfectant should not be disabled. Finally, the computer should be operated in Finder and not Multifinder mode.

Once all the hardware requirements are complete, the tape cassette is inserted into the recorder and the required data found manually. The tape speed should be set at 1.2 cm/s and the tape positioned about 6 seconds before the first peak required. A simple procedure such as counting seconds is usually sufficiently accurate for this purpose. The compiled program is then activated to commence data acquisition.

## *5. The Data Acquisition Program*

A listing of the QuickBasic program used to acquire data from one channel is given in Appendix I. Lines 18 to 24 handle the initialisation of the IEEE interface. The next section down to line 40 initialises program variables before displaying Window2. Lines 44 to 59 display the input window, buttons and default values set by the program. This window is displayed until the pointer is clicked on the CONTINUE button. The next section down to Window3 changes the program variables to reflect what the user has entered on the input window. Window3 displays a progress report on how many records are left to transfer. Statements beginning 'IeeeWrite' send various commands to the IEEE interface, such as telling the recorder to start. For a detailed description of what each command does the Kyowa manuals should be consulted. Line 117 converts the voltage from the recorder to kiloPascals (kPa) and computes the time elapsed since the first record was acquired. The factor 6/5 is used when a backup tape is being scanned. Backup tapes are recorded with an output range of 5 Volts and an input range of 6 Volts, hence the 6/5 conversion. If the original tape is used (although this is not recommended) this factor should be removed and the program recompiled. Lines 129 onwards display a termination window allowing the operator to continue with another record or to quit. The time taken for the previous transfer is also displayed. The subroutines from line 152 onwards are supplied by Iotech and are used to access the IEEE interface.

Appendix II is a listing of a similar program to acquire data when a trigger channel has been recorded as well.



## 6. Operating the Acquisition Program

The Microsoft QuickBasic environment allows compilation of programs to produce an executable Macintosh module. It is this compiled version, executed by double clicking on its Icon, which displays on the screen a window similar to Figure 2. The contents of the five rectangular boxes can be edited using the click and drag function of the Macintosh. The OUTPUT FILE NAME can be any legal Macintosh name. The Press (kPa/volt) box is the calibration value of the pressure transducer used. VTR Range is the voltage input range used when the data were originally acquired, usually somewhere between 0.4 and 6 volts. The No.of Points can be anything from about 2000 to 12 000 but 6000 is a good working compromise between taking a long time and risking missing the peak. A limited number of options is given for the Time interval, this is the interval that the A-to-D Converter will use to acquire data from the recorder. The user can enter 0.5, 1, 2 or 5 ms. These give sample rates of 2000, 1000, 500 and 200 Hz respectively. Any figure greater than 5 is set to 5, any figure less than 0.5 is set to 0.5. Other numbers will be converted to one of the four possible intervals. When converted from the recording speed of 76 cm/s to the playback speed of 1.2 cm/s these intervals convert to real time intervals from 7.9 to 79  $\mu$ s. It is possible to use other values with a small program change but these have been found to be of limited use.

File Edit 16:13:23

PROGRAM PARAMETERS

Alter or Verify these parameters

Click on 'CONTINUE' when ready


OUTPUT FILE NAME	File01.01
Press (kPa/Volt)	7.3
VTR Range	1.4
No.of Points (<= 12000)	6000
Time interval (ms)	.5

CONTINUE

Figure 2: Data input window for data acquisition program.

Clicking on the CONTINUE button on the screen will resume program execution. The recorder should start and after a delay of about 5 seconds the A-D light on the A-to-D converter will come on. The Macintosh screen will change to a window similar to Figure 3. The operator should observe the monitored value in the top right hand corner of the recorder screen to verify that a peak has arrived. If this fails to happen, the operator should wait until the recorder stops and returns to LOCAL mode and should then stop the program by dragging the File menu item on the Macintosh to QUIT. The Macintosh will return to the Finder. The tape should then be repositioned and the program restarted. If the operator is satisfied that the peak has been acquired he can leave the machines for a few minutes until the Macintosh displays a window similar to Figure 4. Clicking on the NEW RECORD button will return the program to the input window. Clicking on the QUIT button will terminate the program.

When a few records have been transferred (for instance, all the channels for one shot) the program can be stopped and the acquired files can be copied to another computer, either via disk or local network. These files can be manipulated while more data is being transferred from the recorder to the computer. It takes 3 or 4 minutes to transfer 6000 points into one file and this time can be used to prepare the already acquired data for graphical display and printing.



*Figure 3: Macintosh window while data is being transferred.*




Figure 4: Termination window for data acquisition program.

## 7. Manipulating the Data

The data already acquired can be quickly scanned by using Igor. It is not within the scope of this Technical Note to describe in detail the use of Igor and only brief explanations are given. On opening Igor the Macintosh screen appears as shown in Figure 5. Pulling down the FILE menu to LOAD WAVES ... LOAD GENERAL TEXT will prompt the user for a filename. On opening the required file the user will be prompted to supply two or three Wave names. Data are stored in the order Times, Pressures then Trigger, if recorded. Pulling down the WAVES menu item to DISPLAY WAVES will prompt the user to select the X and Y variables to plot. The plot should be similar to Figure 6. Editing the data can be accomplished by pulling down the WAVES menu to EDIT WAVES and selecting all the waves displayed in the list. A portion of the data can be kept by highlighting the range in the table that is not wanted and using the CUT function (Command-X on the keyboard) to remove those values. In addition, if times of arrival are required, the arrival time of the trigger can be subtracted from the time values using a simple arithmetic expression. If the trigger time is 31.4 ms and the time data is labelled Times, typing in the Igor Command Window  $\text{Times} = \text{Times} - 31.4$  will alter all the time values by subtracting 31.4. This edited data should then be stored in a Macintosh file in text format. Dragging the FILE menu to SAVE WAVES ... SAVE IGOR TEXT will prompt the user to select the waves to be stored. Only the Time and Pressure should be stored as the trigger is no

longer required. The user will be prompted for a file name to store the data. The Igor program inserts text and unwanted TABS in the file. The text can simply be deleted and the TABS converted to spaces using the Macintosh Desk Accessory McSink. TABS are not handled correctly by Microsoft QuickBasic when reading input files. While editing using McSink, it is also necessary to insert the number of data values as the first line of the data. This is the same as No.of Points in Section 6 and is required by the next program Zerocorr, which is used to compute a mean zero value for the pressure data and operates by prompting the user for up to eight input and output files. The input window is similar to Figure 7. All the rectangular boxes can be edited using Macintosh click and drag procedures. Where less than eight files are available the first input box after the last file should be set to blank. Only files before the blank will be processed. The data is then ready for final output, using Igor to add such things as titles, peak heights and other useful information and to send a copy to a printer. It has been common practice to store the raw data and the Igor graph on floppy disk at this time and also on a backup medium. A 10 Mbyte Jasmine drive has been used for this purpose, however, the Macintosh must be shut down to connect the Jasmine drive or another computer used.




Figure 5: Igor window when first opened.




Figure 6: Igor graph of pressure vs time.

The figure shows a screenshot of a software window titled "PROGRAM PARAMETERS - NO TRIGGER". The window has a menu bar with "File" and "Edit", and a clock in the top right corner showing "16:17:26". Below the title bar, a button labeled "ENTER OR VERIFY FILENAMES, THEN CLICK CONTINUE" is present. The main area contains two columns of file name inputs. The left column is titled "INPUT FILE NAMES" and contains eight text boxes with the following text: "File01.input", "File02.input", "File02.input", "File02.input", "File05.input", "File06.input", "File07.input", and "File08.input". The right column is titled "OUTPUT FILE NAMES" and contains eight text boxes with the following text: "File01.output", "File02.output", "File03.output", "File04.output", "File05.output", "File06.output", "File07.output", and "File08.output". At the bottom of the window, there are two buttons: "QUIT" and "CONTINUE".

Figure 7: Input window for zero correction program.

## ***8. Conclusion***

A simple data acquisition system has been developed from existing equipment with the exception of the Mac/CSI488 interface which was acquired at a small, once only cost. The system is simple to use and, with a copy of this report to hand, personnel should be able to operate it in a relatively short time.

## Appendix I - Listing of Program for One Channel Acquisition

```

00001 '-----
00002 '          DATA TRANSFER FROM VTR TO MACINTOSH - ONE CHANNEL
00003 '-----
00004 '
00005 '    TRANSFER PRESSURE DATA FROM KYOWA VTR TO MACINTOSH
00006 '    STORES DATA POINTS IN A FILE.  THIS VERSION CAN BE COMPILED
00007 '    DATE:          8/5/92
00008 '    FILENAME:      General.1chan
00009 '    AUTHOR:        D.R.KIRK
00010 '    LANGUAGE:      MICROSOFT QUICKBASIC
00011 '    INSTRUMENTS:   KYOWA A-D CONVERTER ADC-116C (IEEE ADDRESS 10)
00012 '                  KYOWA VTR RTP650A (IEEE ADDRESS 15)
00013 '                  MACINTOSH COMPUTER, ANY WITH SCSI PORT
00014 '    INTERFACE:     IOTECH MAC SCSI488 INTERFACE
00015 '-----
00016 '    Ieee Initialization
00017 '    OPTION BASE 0
00018 '    LIBRARY "Device.LIB"          'Library of Basic routines
00019 '    IEEE%=0                      'Interface buffers
00020 '    BufLen%=64
00021 '    Buf$=SPACES(BufLen%)
00022 '    DIM ControlParams$(10)
00023 '-----
00024 '    IeeeInit
00025 '    'Clear and localise AtoD & VTR
00026 '    IeeeWrite "CLEAR 10,15":      IeeeWrite "LOCAL 10,15"
00027 '-----
00028 '-----
00029 '          'Restart here if more data to analyse
00030 '    Newstart:
00031 '          'Set initial values
00032 '    Npts%=6000                    'Number of data values
00033 '    Outrange=5!                  'VTR output range (Volts)
00034 '    Playback=1.2                 'VTR speed cm/s
00035 '    Record=76!                  'VTR speed cm/s
00036 '    Adinterval=.5               'AD Converter interval (ms)
00037 '    Range=5!                    'AD Voltage range
00038 '    channel$="File01.01"        'File name for data storage
00039 '    chncal=7.3                  'Pressure kPa/Volt
00040 '    Inrange=1.4                 'VTR input range
00041 '-----
00042 '    Window2:                    ' display initial values in edit fields
00043 '-----
00044 '    WINDOW 2 , "PROGRAM PARAMETERS", (0,40)-(525,336)
00045 '    BUTTON 5,1, "Alter or Verify these parameters", (140,20)-(375,35),1
00046 '    BUTTON 6,1, "Click on 'CONTINUE' when ready", (140,40)-(375,55),1
00047 '    EDIT FIELD 11, "Press (kPa/Volt)", (105,120)-(260,135),3
00048 '    EDIT FIELD 12, STR$(chncal), (270,120)-(400,135)
00049 '    EDIT FIELD 13, "VTR Range", (105,140)-(260,155),3
00050 '    EDIT FIELD 14, STR$(Inrange), (270,140)-(400,155)
00051 '    EDIT FIELD 15, "No. of Points (<= 12000)", (105,160)-(260,175),3
00052 '    EDIT FIELD 16, STR$(Npts%), (270,160)-(400,175)
00053 '    EDIT FIELD 17, "Time interval (ms)", (105,180)-(260,195),3
00054 '    EDIT FIELD 18, STR$(Adinterval), (270,180)-(400,195)
00055 '    EDIT FIELD 9, "OUTPUT FILE NAME", (105,100)-(260,115),3
00056 '    EDIT FIELD 10, channel$, (270,100)-(400,115)
00057 '    BUTTON 1,1, "CONTINUE", (300,245)-(370,260),1
00058 '    ACTIVITY=DIALOG(0)
00059 '    WHILE ACTIVITY<>1:ACTIVITY=DIALOG(0):WEND
00060 '-----
00061 '    Extract altered values from edit fields
00062 '-----
00063 '    Inrget%=INT(Range)/4
00064 '    channel$=EDIT$(10)
00065 '    chncal=VAL(EDIT$(12))
00066 '    Inrange=VAL(EDIT$(14))
00067 '    Npts%=VAL(EDIT$(16))
00068 '    Adinterval=VAL(EDIT$(18))
00069 '    Runtime=Npts%*Adinterval*.001+2
00070 '    Timeint=Adinterval*Playback*.001/Record

```

```

00071 IF Adinterval<1 THEN
00072     Clk$="CLK6"
00073     GOTO Window3
00074 END IF
00075 IF Adinterval<1.8 THEN
00076     Clk$="CLK7"
00077     GOTO Window3
00078 END IF
00079 IF Adinterval<4 THEN
00080     Clk$="CLK8"
00081     GOTO Window3
00082 END IF
00083 Clk$="CLK9"
00084 LOCATE 15,10
00085 Window3: '-----
00086 '                Transfer data from VTR to A to D Converter
00087
00088 Tzt=TIMER
00089 WINDOW CLOSE 2
00090 WINDOW 3,*** TRANSFERRING DATA ***,(0,40)-(525,336)
00091 BUTTON 1,1,"DATA FROM VTR TO AD CONV", (170,25)-(355,40),1
00092 IeeeWrite "OUTPUT 15;PF,EX" '15 is GPIB address of VTR, start VTR
00093 '     Reset AtoD, storage mode, rate
00094 IeeeWrite "OUTPUT 10;RST,MOD1,"+Clk$ '10 is address of AtoD converter
00095 '     Set for channel 1, range 1 or 5 volts
00096 IeeeWrite "OUTPUT 10;CHN001001,RNG"+MID$(STR$(Irange%),2,1)
00097 ZTIME=TIMER
00098 WHILE TIMER-ZTIME<5: WEND 'Delay A to D conversion to allow VTR
00099 IeeeWrite "OUTPUT 10;TRM0,STA" 'to come up to speed
00100 ZTIME=TIMER
00101 WHILE TIMER-ZTIME<Runtime: WEND
00102 IeeeWrite "OUTPUT 10;STP" 'Stop AtoD conversion
00103 IeeeWrite "OUTPUT 15;SP,EX" 'Stop VTR
00104 IeeeWrite "CLEAR15": IeeeWrite "LOCAL15"
00105 OPEN channel$ FOR OUTPUT AS #1
00106 BUTTON 1,1,"DATA FROM AD CONV TO FILE", (170,25)-(355,40),1
00107 '     Set AtoD to start reading at chn1, all channels, number of points
00108 IeeeWrite "OUTPUT 10;SPN1,RCH0,RDD"+MID$(STR$(Npts%),2)
00109 LOCATE 6,23
00110 PRINT " Records remaining"
00111 limit%=100
00112 k%=0
00113     WHILE k%<Npts%
00114         FOR i%=1 TO limit%
00115             k%=k%+1
00116             A$="": IeeeWrite "ENTER 10#6": IeeeRead A$: A1%=VAL(A$)
00117             Pressure=A1%*chncl*Range*Irange/(2000*Outrange)
00118             Time1=(k%-1)*Timeint
00119             PRINT #1,Time1,Pressure*6/5 '6/5 only used for backup tape
00120         NEXT i%
00121         LOCATE 6,40
00122         PRINT Npts%-k%; " "
00123         WEND
00124 IeeeWrite "CLEAR10": IeeeWrite "LOCAL10"
00125 CLOSE #1
00126 LOCATE 6,40
00127 PRINT " "
00128 '-----
00129 Finis:
00130 TOTIME=TIMER-Tzt
00131 BUTTON 3,1,"WARNING!, REWIND VTR TO START POINT BEFORE
        PROCEEDING", (50,230)-(450,245),1
00132 BUTTON 1,1,"QUIT", (80,260)-(120,275),1
00133 BUTTON 2,1,"NEW RECORD", (150,260)-(250,275),1
00134 LOCATE 12,20
00135 PRINT "TIME TAKEN = ";TOTIME;" SECONDS"
00136 ACTIVITY=DIALOG(0)
00137     WHILE ACTIVITY<>1
00138         ACTIVITY=DIALOG(0)
00139     WEND
00140 ON DIALOG(1) GOTO Termin,Newst,Finis
00141 Newst:
00142     BUTTON CLOSE 1
00143     BUTTON CLOSE 2

```



```

00144 BUTTON CLOSE 3
00145 WINDOW CLOSE 3
00146 GOTO Window2
00147 Termin:
00148 WINDOW CLOSE 3
00149 Quits:
00150 END
00151 '-----
Note: The following subroutines appear at the end of both programs, although
      they are omitted from the second listing.
00152 ' Ieee SubRoutines required by Iotech MacSCSI488 bus controller
00153 SUB IeeeInit STATIC
00154     SHARED IEE%
00155     written%=0
00156     OpenDriver ".IEE",IEE%
00157     IeeeRestart
00158     IeeeWrite "RESET"
00159     IeeeWrite "EOL IN NONE"
00160 END SUB
00161
00162 SUB IeeeFinish STATIC
00163     SHARED IEE%
00164     CloseDriver IEE%
00165 END SUB
00166
00167 SUB IeeeRestart STATIC
00168     SHARED IEE%,ControlParams%()
00169     Control IEE%,256,ControlParams%(0)
00170 END SUB
00171
00172 SUB IeeeState(State%) STATIC
00173     SHARED IEE%,ControlParams%()
00174     Status IEE%,0,ControlParams%(0)
00175     State%=ControlParams%(0)
00176 END SUB
00177
00178 SUB IeeeActive(State%) STATIC
00179     SHARED IEE%,ControlParams%()
00180     Status IEE%,0,ControlParams%(0)
00181     State%=ControlParams%(1)
00182 END SUB
00183
00184 SUB IeeeWriteB(text%) STATIC
00185     SHARED IEE%
00186     written%=0
00187     FSWrite IEE%,text$,written%
00188 END SUB
00189
00190 SUB IeeeWrite(text%) STATIC
00191     IeeeWriteB(text$+CHR$(13))
00192 END SUB
00193
00194 SUB IeeeRead(result%) STATIC
00195     SHARED IEE%,Buf$
00196     ReadIn%=0
00197     FSRead IEE%,Buf$,ReadIn%
00198     result$=LEFT$(Buf$,ReadIn%)
00199     IeeeState State%
00200     WHILE State%<>0
00201         FSRead IEE%,Buf$,ReadIn%
00202         result$=result$+LEFT$(Buf$,ReadIn%)
00203         IeeeState State%
00204     WEND
00205 END SUB

```

## Appendix II - Listing of Program for Two Channel Acquisition

```

00001 '-----
00002 '          DATA TRANSFER FROM VTR TO MACINTOSH - TWO CHANNEL
00003 '-----
00004 '
00005 '   TRANSFER KYOWA VTR DATA AND TRIGGER DATA TO A MACINTOSH FILE
00006 '       DATE:          8/5/92/91
00007 '       FILENAME:      General.trig
00008 '       AUTHOR:        D.R.KIRK
00009 '       LANGUAGE:      MICROSOFT QUICKBASIC
00010 '       INSTRUMENTS:   KYOWA A-D CONVERTER ADC-116C
00011 '                       KYOWA VTR RTP650A
00012 '                       MACINTOSH SE OR SE/30
00013 '       INTERFACE:     MACSCSI 488 - IOTECH
00014 '-----
00015 '       Ieee Initialization
00016 '   OPTION BASE 0
00017 '   LIBRARY "Device.LIB"
00018 '   IEEE%=0
00019 '   BufLen%=64
00020 '   Buf$=SPACES$(BufLen%)
00021 '   DIM ControlParams$(10)
00022 '-----
00023 '   IeeeInit
00024 '   IeeeWrite "CLEAR 10,15":          IeeeWrite "LOCAL 10,15"
00025 '                                   'Clear          and localise AtoD & VTR
00026 '-----
00027 '-----
00028 '   channels$="File01.01"          'File name for data storage
00029 '   Npts%=6000
00030 '   Inrange=6          'VTR input range
00031 '   Outrange=5!        'VTR output range
00032 '   Playback=1.2       'VTR cm/s
00033 '   Record=76!         'VTR cm/s
00034 '   Adinterval=1!      'AD Converter (ms)
00035 '   Range=5!           'AD Voltage range
00036 '   Irnge%=INT(Range)/4
00037 '   chncal=7!
00038 '   Newstart:          'Restart here if more data to analyse
00039 '   Window2:
00040 '   WINDOW 2,"PROGRAM PARAMETERS",(0,40)-(525,336)
00041 '   BUTTON 5,1,"Alter or Verify these parameters",(140,20)-(375,35),1
00042 '   BUTTON 6,1,"Click on 'CONTINUE' when ready",(140,40)-(375,55),1
00043 '   EDIT FIELD 11,"Pressure (kPa/Volt)",(105,120)-(250,135),3
00044 '   EDIT FIELD 12,STR$(chncal),(260,120)-(400,135)
00045 '   EDIT FIELD 13,"VTR Range",(105,140)-(250,155),3
00046 '   EDIT FIELD 14,STR$(Inrange),(260,140)-(400,155)
00047 '   EDIT FIELD 15,"Time Interval (ms)",(105,160)-(250,175),3
00048 '   EDIT FIELD 16,STR$(Adinterval),(260,160)-(400,175)
00049 '   EDIT FIELD 17,"No.of Points (<8000)",(105,180)-(250,195),3
00050 '   EDIT FIELD 18,STR$(Npts%),(260,180)-(400,195)
00051 '   EDIT FIELD 9,"OUTPUT FILE NAME",(105,100)-(250,115),3
00052 '   EDIT FIELD 10,channels$,(260,100)-(400,115)
00053 '   BUTTON 1,1,"CONTINUE",(300,230)-(370,245),1
00054 '   ACTIVITY=DIALOG(0)
00055 '   WHILE ACTIVITY<>1:ACTIVITY=DIALOG(0):WEND
00056 '   channels$=EDIT$(10)
00057 '   chncal=VAL(EDIT$(12))
00058 '   Inrange=VAL(EDIT$(14))
00059 '   Adinterval=VAL(EDIT$(16))
00060 '   Npts%=VAL(EDIT$(18))
00061 '   Runtime=Npts%*Adinterval*.001+2
00062 '   Timeint=Adinterval*Playback/Record
00063 '   IF Adinterval<1 THEN
00064 '       Clk$="CLK6"          ' Ad interval = 0.5 ms Real time 7.9 us
00065 '       GOTO Window3
00066 '   END IF
00067 '   IF Adinterval<1.8 THEN
00068 '       Clk$="CLK7"          ' Ad interval = 1.0 ms Real time 15.8 us
00069 '       GOTO Window3

```

```

00070 END IF
00071 IF Adinterval<4 THEN
00072   Clk$="CLK8" ' Ad interval = 2.0 ms Real time 31.7 us
00073   GOTO Window3
00074 END IF
00075 Clk$="CLK9" ' Ad interval = 5.0 ms Real time 78.9 us
00076 LOCATE 15,10
00077 Window3: '-----
00078 Tzt=TIMER
00079 WINDOW CLOSE 2
00080 WINDOW 3,"*** TRANSFERRING DATA ***", (0,40)-(525,336)
00081 BUTTON 1,1,"DATA FROM VTR TO AD CONV", (170,25)-(355,40),1
00082 IeeeWrite "OUTPUT 15;PF,EX"
00083 IeeeWrite "OUTPUT 10;RST,MOD1,"+Clk$ ' Storage Mode,
Interval
00084 IeeeWrite "OUTPUT 10;CHN001002,RNG"+MID$(STR$(Irnge%),2,1)
00085 ztime=TIMER
00086 WHILE TIMER-ztime<5: WEND 'delay A to D start by 5
second
00087 IeeeWrite "OUTPUT 10;TRM0,STA" 'Start A to D
conversion
00088 ztime=TIMER
00089 WHILE TIMER-ztime<Runtime: WEND
00090 IeeeWrite "OUTPUT 15;SP,EX" 'Stop VTR
00091 IeeeWrite "OUTPUT 10;STP" 'Stop A to D
converter
00092 IeeeWrite "CLEAR15": IeeeWrite "LOCAL15" 'Clear and localise VTR
00093 BUTTON 1,1,"DATA FROM AD CONV TO FILE", (170,25)-(355,40),1
00094 IeeeWrite "OUTPUT 10;SPN1,RCH0,RDD"+MID$(STR$(Npts%),2)
' Set A to D for data output
00095 OPEN channel$ FOR OUTPUT AS #1
00096 LOCATE 4,23
00097 PRINT " Records remaining"
00098 limit%=100
00099 k%=0
00100 WHILE k%<Npts%
00101   FOR i%=1 TO limit%
00102     k%=k%+1
00103     A$="": IeeeWrite "ENTER 10#6": IeeeRead A$:
A1%=VAL(A$)
00104     B$="": IeeeWrite "ENTER 10#6": IeeeRead B$:
B1%=VAL(B$)
00105     Pressure=A1%*chncal*Range*Inrange/(2000*Outrange)
00106     PRINT #1,k%*Timeint*.001;CHR$(9);Pressure*6/5;CHR$(9);B1%
00107     '6/5 only for backup
tape
00108     NEXT i%
00109     LOCATE 4,40
00110     PRINT Npts%-k%; " "
00111   WEND
00112   LOCATE 4,40
00113   PRINT " "
00114   IeeeWrite "CLEAR10": IeeeWrite "LOCAL10" 'Clear and localise
AtoD
00115   CLOSE #1
00116   '-----
00117   Finis:
00118   TOTIME=TIMER-Tzt
00119   BUTTON 6,1,"WARNING!, REWIND VTR TO START POINT BEFORE
PROCEEDING", (50,230)-(450,245),1
00120   BUTTON 1,1,"QUIT", (80,260)-(120,275),1
00121   BUTTON 2,1,"NEW RECORD", (150,260)-(250,275),1
00122   LOCATE 12,20
00123   PRINT "TIME TAKEN = ";TOTIME;" SECONDS"
00124   ACTIVITY=DIALOG(0): WHILE ACTIVITY<>1: ACTIVITY=DIALOG(0): WEND
00125   ON DIALOG(1) GOTO Termin,Newst
00126   Newst:
00127   BUTTON CLOSE 1
00128   BUTTON CLOSE 2
00129   BUTTON CLOSE 6
00130   WINDOW CLOSE 3
00131   GOTO Newstart
00132   Termin:
00133   WINDOW CLOSE 3

```

00134 Quits:  
00135 END  
00136 -----  
00137 ' Ieee functions as for one channel listing

REPORT NO.  
MRL-TN-628AR NO.  
AR-008-234REPORT SECURITY CLASSIFICATION  
Unclassified

## TITLE

Data acquisition from blast overpressure trials

AUTHOR(S)  
D.R. KirkCORPORATE AUTHOR  
DSTO Materials Research Laboratory  
PO Box 50  
Ascot Vale Victoria 3032REPORT DATE  
March, 1993TASK NO.  
ARM 89/178SPONSOR  
ArmyFILE NO.  
G6/4/8-4332REFERENCES  
-PAGES  
21

CLASSIFICATION/LIMITATION REVIEW DATE

CLASSIFICATION/RELEASE AUTHORITY  
Chief, Explosives Ordnance Division

## SECONDARY DISTRIBUTION

Approved for public release

## ANNOUNCEMENT

Announcement of this report is unlimited

## KEYWORDS

User manual

Data processing

Blast overpressure

## ABSTRACT

A Macintosh computer has been used to acquire data from blast overpressure trials on various weapons. The computer is connected to a multiple channel FM data recorder via a MacSCS1488 bus controller, allowing the computer to control the recorder and to acquire data from it through an analog to digital converter. Detailed instructions are given for connecting the hardware and operating the software involved.

Data Acquisition from Blast Overpressure Trials

D.R. Kirk

(MRL-TN-628)

DISTRIBUTION LIST

Director, MRL  
Chief, Explosives Ordnance Division  
Mr L. Sheppard  
Mr D.R. Kirk  
MRL Information Service

Chief Defence Scientist (for CDS, FASSP, ASSCM) (1 copy only)  
Director, Surveillance Research Laboratory  
Director (for Library), Aeronautical Research Laboratory  
Director, Electronics Research Laboratory  
Head, Information Centre, Defence Intelligence Organisation  
OIC Technical Reports Centre, Defence Central Library  
Officer in Charge, Document Exchange Centre (8 copies)  
Army Scientific Adviser, Russell Offices  
Air Force Scientific Adviser, Russell Offices  
Navy Scientific Adviser, Russell Offices - data sheet only  
Scientific Adviser, Defence Central  
Director-General Force Development (Land)  
Senior Librarian, Main Library DSTOS  
Librarian, MRL Sydney  
Librarian, H Block  
UK/USA/CAN ABCA Armies Standardisation Rep. c/- DGAT (8 copies)  
Librarian, Australian Defence Force Academy  
Counsellor, Defence Science, Embassy of Australia - data sheet only  
Counsellor, Defence Science, Australian High Commission - data sheet only  
Scientific Adviser to DSTC, C/- Defence Adviser - data sheet only  
Scientific Adviser to MRDC, C/- Defence Attache - data sheet only  
Head of Staff, British Defence Research and Supply Staff (Australia)  
NASA Senior Scientific Representative in Australia  
INSPEC: Acquisitions Section Institution of Electrical Engineers  
Head Librarian, Australian Nuclear Science and Technology Organisation  
Senior Librarian, Hargrave Library, Monash University  
Library - Exchange Desk, National Institute of Standards and Technology, US  
Exchange Section, British Library Document Supply Centre  
Periodicals Recording Section, Science Reference and Information Service, UK  
Library, Chemical Abstracts Reference Service  
Engineering Societies Library, US  
Documents Librarian, The Center for Research Libraries, US

Dr R. Glass, EOD-S  
Dr H.T. French, EOD-S  
Mr J.K. Leach, EOD-S  
Mr P.L. Fagan, EOD-S  
Mr K.L. Krishnamoorthy, EOD-S  
Mr G. Teague, EOD-S  
Dr S. Stephenson, EOD-S  
Mr J.S. Howe, SSD Melbourne